

# SAFE POWERFUL SPEEDY

## UVC Sterilization Box GUV-C105



### Excellent performance


With high-output and powerful ultraviolet LED (UVC: wavelength 260-280nm), our high-efficiency sterilizer kills 99.9% of the bacteria and viruses on your belongings in just 90 seconds.

### Chemical-free

Not to worry about chemical residues as it doesn't require the use of any chemical substances such as alcohol, etc. Moreover, the UV lamps turn off as soon as you open the lid of the box to ensure your safety and security

### 360° Sterilization

Offers wide storage and evenly sterilizes all the surface with 360° LED radiation



**GaN Power Technology Co., Ltd.**

## Use:

### In your home and office

- ✓ Keep it at entrance or in your office
- ✓ Sterilize your belongings whenever you reach home or periodically

### At your business

- ✓ Display eyewear, jewelry, watches after trial
- ✓ Combs & scissors after every use at beauty salons
- ✓ Locker keys at sports complex/gymnasiums
- ✓ Rental audio guide at museums

Test Report

萬能科技大學生物技術系  
Biotechnology Department, Vanung University  
32045 No.1 Van-Nung Rd., Chung-Li, Tao-Yuan, Taiwan, R.O.C.  
TEL: 886-03-4515811, ext 340

分析檢驗報告  
TEST REPORT

委託單位：鴻輝科技股份有限公司 VI2020059901A  
Name of Address Client: SMS  
Sample designation (標品名稱)：功能性 LED (GPT 二代)  
Test Item (試驗項目)：Filtration activity Analysis (過濾檢測)  
Date of Acceptance (收件日期)：05/08/2020  
Date of Issue (報告日期)：05/09/2020

測試項目	Bacteria 菌株	Results 結果			
		接種量 cfu/ml	菌落直徑 mm	過濾率 Activity	照射秒數
<input checked="" type="checkbox"/>	<i>Klebsiella pneumoniae</i> (肺炎桿菌 ATCC 4962)	3.5x10 <sup>7</sup>	200	>99.9%	48 (90min)

Note: 1. 檢驗方法請參考 American Public Health Association, Section AATCC  
2. 對於報告內容請參考 AATCC30, AATCC147

※注意事項 Notes:  
1. 本報告僅針對此次送檢之樣品為分析檢驗結果，僅供委託者內部參考之用，不得從事廣告、公開、商業等行為或作其他用途之使用。  
This report only accounts for the sample(s) sent and should not be used for any commercial advertisement and public statement.  
2. 本報告將禁止以任何方式翻印或複製且不得偽造或篡改。  
This report is invalid without being stamped and dated, or if obliterated.  
3. 如有本報告內容任何問題，歡迎隨時向本單位聯絡查詢。  
If any question about this report, please contact with us.

Lab Leader: D. Grey 實驗室主任:

Verified results by Vanung University (Taiwan)  
Eliminated more than 99.9% Klebsiella  
Pneumoniae bacteria in just 48 seconds

## UVC Sterilization box specifications

Detail	Specification
Model	GUV-C105
Sterilization time	90 seconds
UV wavelength	260-280nm
UV output	Total 105mW (upper lights 100mW, lower lights 5mW)
Interior size	20.0x9.0x11.0 cm
Power input	5VDC 2A
Power interface	USB Type-C
Accessories	User manual, AC adapter, Storage cum Carry bag
Exterior size and weight	22.8x13.9x13.8 cm 965g
Product warranty	1 year
Made in	Taiwan

\*Please read the user instruction manual carefully before use

Manufactured by GaN Power Technology Co., Ltd.

Address: 7F.-10, No. 262, Sec. 2, Henan Rd., Xitun Dist.,  
Taichung City 40744, Taiwan (R.O.C.)

TEL: +886-4-24528666

E-mail: service@ganpowergroup.com

Sales partner:

2007A